

# Kindera Replaceable Part - Caster

## Install Instructions

### Tools Required

- Wrench
- Rubber Mallet

### Installation

1. Lay the recliner on its back.
2. Locate the caster stem, washer, and caster (Figure A).
3. Insert the washer on the threaded side of the caster stem.
4. Thread the caster stem into the hole (Detail A, Figure A).
5. Using a wrench, tighten the caster stem into place.
6. Line up the caster with the caster stem and use a rubber mallet to seat onto the stem (Detail B and Detail C, Figure B).
7. The caster is locked into place when the caster is pulled on and does not come off.
8. The caster should be able to spin 360 degrees freely.
9. Flip the recliner up right (Figure C).

Figure A




Figure B




Figure C

